

It's a
feel-good
thing



G·U·M[®]

**GUIDA PRODOTTI
FARMACIA**

2025



INTERDENTALI

SOFT-PICKS[®] PRO



689M30 GUM SOFT-PICKS[®] PRO S
987318072



690M12 GUM SOFT-PICKS[®] PRO M
987318096



690M30 GUM SOFT-PICKS[®] PRO M
987318084



691M30 GUM SOFT-PICKS[®] PRO L
987318108



BI-DIRECTION



2014 GUM BI-DIRECTION ISO 0
983696497



2114 GUM BI-DIRECTION ISO 1
978503896



2314 GUM BI-DIRECTION ISO 2
978503908



2614 GUM BI-DIRECTION ISO 3
978503910



2714 GUM BI-DIRECTION ISO 4
982543985



SOFT-PICKS[®] MINTY COOL MINT FLAVOUR



670BT40M1 GUM SOFT-PICKS[®]
MINTY M 40 pcs
978503860



670BT80M1 GUM SOFT-PICKS[®]
MINTY M 80 pcs
985824059



669BT40M1 GUM SOFT-PICKS[®]
MINTY S 40 pcs
985824010



671BT40M1 GUM SOFT-PICKS[®]
MINTY L 40 pcs
985824034



TRAV-LER



1312 0,6 mm
935236253



1314 0,8 mm
935236265



1412 0,9 mm
935236277



1414 1,1 mm
935236289



1512 1,2 mm
972036305



1514 1,3 mm
972036317



1612 1,4 mm
972036329



1614 1,6 mm
935236303



1618 2,0 mm
972036331



1619 2,6 mm
980806956





SPAZZOLINI PER ADULTI



NEW



Prevenzione giornaliera

CLASSIC

409 GUM Classic
compatto morbido
902220298



410 GUM Classic
compatto medio
932121458



411 GUM Classic
regolare morbido
902220312



Prevenzione giornaliera

EASY CLEAN



135 GUM EASY-CLEAN
compatto Medio
989332642



GUM 305



305 GUM 305
compatto duro
906664545



Ipersensibilità

Sonic SENSITIVE



4101 GUM SONIC
SENSITIVE
982938615



4111 GUM SONIC
SENSITIVE refills (2x)
982938627



PRO SENSITIVE



510 GUM PRO
SENSITIVE
986873533



Prevenzione giornaliera avanzata

Sonic DAILY



4100MWH2 GUM
SONIC DAILY WHITE
982938538



4110MWH2 GUM
SONIC DAILY
WHITE REFILLS (2x)
982938603



4100MBK2 GUM
SONIC DAILY BLACK
982938490



4110MBK2 GUM
SONIC DAILY BLACK
refills (2x)
982938591



PRO



525
GUM PRO
compatto
morbido
934323724



526
GUM PRO
regolare
medio
934323763



528
GUM PRO
compatto
medio
934323775





SPAZZOLINI SPECIALI E ACCESSORI

Prevenzione giornaliera

Classic 311



311 GUM Classic
slender morbido
902220324



ActiVital[®]



583 GUM Activital[®]
medio
937463141



585 GUM Activital[®]
morbido
930007327



Post-operatorio



317 GUM
Post-
Operation
903026793



Denti bianchi
Original WHITE[®]



563 GUM
Original
White
927383873



Ortodonzia
ORTHO[®]



124 GUM Ortho
ortodontico
902222328



Protesi dentali



201 GUM Denture
Brush per protesi
902222850



Puliscilingua



760 GUM
HaliControl[®]
puliscilingua
903760991



Monociuffo



308 GUM End-Tuft
monociuffo conico
902222847



Kit da viaggio



156 GUM
Travel Kit
902311606



Pastiglie rivela placca



824 GUM
Red-Cote[®]
pastiglie
rivela placca
902224195





FILI



3500
GUM TWISTED FLOSS
984881019



7 630019 805107



2030M
GUM Expanding Floss
900415100



0 070942 302401



3220M
GUM ORTHO Floss
933015238



7 630019 900935



2000
GUM Easy-Floss
935689808



0 070942 020008



2040
GUM Original White
904802749



0 070942 303200



1155
GUM Butlerweave[®] cerato
906668076



0 070942 011556



890
GUM Easy Flossers NEW
980501629



7 630019 903783



PROTEZIONE GENGIVE AZIONE MANTENIMENTO

PAROEX[®]



1702 GUM PAROEX[®]
0,06% CHX + 0,05%
CPC collutorio 500 ml
938123142



526 GUM PRO regolare
medio 934323763



1750 GUM PAROEX[®]
0,06% CHX + 0,05%
CPC dentifricio 75 ml
938123167



690M30 GUM
SOFT-PICKS[®] PRO
987318084



PROTEZIONE GENGIVE AZIONE POTENZIATA

PAROEX[®]



1784 GUM PAROEX[®] 0,12% CHX
+ 0,05% CPC collutorio 300 ml
907048247



526 GUM PRO
regolare medio
934323763



1790 GUM PAROEX[®] 0,12% CHX
+ 0,05% CPC dentifricio 75 ml
908008927



1612 GUM
TRAV-LER[®] 1,4 mm
972036329



1785 GUM PAROEX[®] 0,20% CHX
collutorio 300 ml
931507394



DENTI SENSIBILI

SensiVital⁺



6081 GUM SensiVital⁺
collutorio 300 ml
975508971



510 GUM PRO Sensitive
spazzolino 986873533



6070 GUM SensiVital⁺
dentifricio 75 ml
975508957



4101 P GUM SONIC
SENSITIVE 982938615



IGIENE ORALE QUOTIDIANA

ActiVital[®]



6060 GUM ActiVital[®]
collutorio 300 ml
971347075



4100MWH2 GUM
SONIC DAILY
975508983



6050 GUM ActiVital[®]
dentifricio 75 ml
971347063



670BT40M1 GUM
SOFT-PICKS[®] MINTY M
978503860





IGIENE ORALE BIMBI

BABY (0-2)



213 GUM Baby Monster (0-2 anni)
930007416



0 070942 1255 12

KIDS (3-6)



3000 GUM Kids Monster (3+ anni)
dentifricio 50 ml
930007354



0 070942 304153

901 GUM Kids Monster (3-6 anni)
spazzolino
971347087



7 630019 902557

JUNIOR (6+)



SI ILLUMINA

3022 GUM Junior (6+ anni) collutorio
300 ml
930007404



7 630019 902892

3004 GUM Junior Monster (6+ anni)
dentifricio 50 ml
930007392



0 070942 304160

902 GUM Junior Monster (7-9 anni)
spazzolino
971347099



7 630019 902564

903 GUM Light UP Monster spazzolino
971347101



7 630019 902595



BOCCA SECCA

HYDRAL'



6000 GUM HYDRAL' gel
50 ml
934625878



7 630019 901727

6010 GUM HYDRAL' spray
50 ml
934625854



7 630019 901758

6020 GUM HYDRAL'
dentifricio 75 ml
934625892



7 630019 901741



SBIANCANTE ANTI-MACCHIA

Original WHITE



1745 GUM Original White
dentifricio 904396367



0 070942 303132

563 GUM Original White
spazzolino 927383873



7 630019 902304

2040 GUM Original White
floss 904802749



0 070942 303200



AFTE

Afta Clear



2420 GUM AftaClear' spray 15 ml
970435830



7 630019 902229

2400 GUM AftaClear' gel 10 ml
970435816



7 630019 902205

2410 GUM AftaClear' collutorio 120 ml
970435828



7 630019 902212



IGIENE ORALE BIOLOGICA



7020 GUM BIO dentifricio
biologico 75 ml 981472602



7030 GUM BIO Collutorio biologico
300 ml 988028369



ORTODONZIA

ORTHO®



3090 GUM ORTHO
collutorio 300 ml
931156867



3080 GUM ORTHO
dentifricio
931156879



124 GUM ORTHO
spazzolino
90222328



723 GUM ORTHO cera
902224144



3220M GUM ORTHO
Floss
933015238



PROBIOTICI

PerioBalance®



7010 GUM PerioBalance® compresse
905306775



G·U·M®

Scarica qui il
catalogo dei
prodotti:



SUNSTAR ITALIANA S.R.L.

[f](https://www.facebook.com/sunstargumitalia) /sunstargumitalia

[@](https://www.instagram.com/gum_italia) /gum_italia

[in](https://www.linkedin.com/company/sunstar-italiana-s.r.l.) company/sunstar-italiana-s.r.l.

[▶](https://www.youtube.com/sunstargumitalia) /sunstargumitalia

SO125CATBREVEPH